



progress in surgery

UVDID
HYBRID

Hysteroscopy Set




HYSTEROSCOPY SET

HYSTEROSCOPY

Office Hysteroscopy is getting more and more important.
Diagnosis and treatment at the same time of uterus pathologies in a safe and minimal invasive procedure becomes the new standard.



DIAGNOSTIC HYSTEROSCOPE SET




The 4mm Hysteroscope offers all advantages for intrauterine, atraumatic diagnostic.

HYSTEROSCOPE
Ø 4.0mm, 30°, 300mm working length
351-804-030




DIAGNOSTIC (INNER-)SHEATH
Ø 5.4mm, with stopcock, inclusive obturator
260mm working length
351-500-054




OUTER SHEATH
Ø 6.5mm, with stopcock
240mm working length
351-500-065




MIDI OPERATION-HYSTEROSCOPE SET Ø 2.9mm




HYSTEROSCOPE
Ø 2.9mm, 30°, 300mm working length
351-829-030



OPERATION (INNER-)SHEATH
with stopcock, working channel for 5Charr. instruments
extended working length: 216mm
351-400-143



OUTER SHEATH
Ø 5mm, pear-shaped, with stopcock
extended working length: 216mm
351-400-150



HYBRID

With the RZ Hybrid Resectoscope, we provide the latest technology which allows bipolar or monopolar working procedures with one single instrument set. The change between bipolar or monopolar can be done with minor adaptations. The HYBRID Electrodes can be used bipolar or monopolar. Besides the change of the irrigation fluid Glycine/Sorbitol is used during monopolar application and NaCl solution for bipolar application, just the cable needs to be changed.



HYBRID HYSTEROSCOPE SET

HYSTEROSCOPE

Ø 4.0mm, 12°, 300mm working length

351-804-012

PASSIVE WORKING ELEMENT


for HYBRID application

253-000-400

CONTINUOUS FLOW SHEATH

26 FR., 360° rotatable, oblique ceramic tip
Quick Lock mechanism









253-000-256




HYSTEROSCOPY SET

ELECTRODES

Specially designed for RZ HYBRID Resectoscope

		
	Loop, straight	253-400-200
	Loop, curved	253-400-201
	Loop, curved, strong wire	253-400-202
	Collins Coagulation Electrode, pointed	253-400-203
	Coagulation Electrode, conical	253-400-204
	Coagulation Electrode, ball shaped, Ø 3mm	253-400-205
	Roller-Electrode, Ø 3mm	253-400-211

NOTES





progress in surgery

RZ Medizintechnik GmbH

Unter Hasslen 20
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

