



progress in surgery

STERNUM PLATE & WIRE CUTTER



STERNUM PLATE & WIRE CUTTER

Sternal Implants such as plates, wires or similar devices are used for closing and fixing the sternum after Thoracotomy. In several cases the implants needs to be removed. That could be an elective surgery or even emergency surgery.

For cutting sternal plates or wire we have designed a dedicated Implant cutter. This Sternum Implant cutter with a special cutting tip design and material allows cutting of implants directly on the bone surface.

Even materials thicker than 1.7mm can be cut directly on the sternal bone surface.





STERNUM Implant Cutter

230mm, 9"


215-010-040



STERNUM Implant Cutter coated
— special coated tip

230mm, 9"

215-010-042





progress in surgery

RZ Medizintechnik GmbH

Unter Hasslen 20
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

