




progress in surgery

ArthVANCE[®]

Suture Passer





ARTHVANCE® SUTURE PASSER


Rotator Cuff Repair

Simultaneously one-handed passing and retrieving of suture

- Smooth and easy suture passing
- Efficient tissue bite (wide jaw opening and deep bite)
- Minimal perforation of tissue (small hole)
- Reusable instrument and reliable construction
- Simple assembling/disassembling




Suture Passer with trap door frontal opening
4-finger trigger, **without ratchet**
170mm working length 225-820-217




Suture Passer with trap door frontal opening
4-finger trigger, **with ratchet**
170mm working length 225-820-213



Needle with lateral opening
0.3 x 240mm 225-820-219



Automatic trap loop hold



Slim design



Ratchet mechanism
(225-820-213)



Strong grip
4-finger trigger



progress in surgery

RZ Medizintechnik GmbH

Unter Hasslen 20
78532 Tuttlingen, Germany

Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

